OSMODEMI

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In the OSMODEMI demineralizers, the water to be treated is fed sequentially to an activated charcoal cartridge, a 5 micron safety filter and a RO membrane module. Pure water leaving the RO module (i.e. RO water) is eventually fed to a cartridge packed with mixed-bed ion exchange resins that further increase water purity to $0.5 \mu\text{S/cm}$ (i.e., $> 2 \text{ M}\Omega/\text{cm}$). In the ion exchange cartridge, about 4 liters of high efficiency, long-life mixed resins are packed that ensure long production times. A digital display on the control panel connected to an in-line conductivity meter continuously displays the produced water conductivity and permits control of its purity. Increasing water conductivity signals exhaustion of the resin salt removal capacity. Cartridges are very easily replaced.

The OSMODEMI demineralizers are supplied in two configurations. Demineralizers in the STANDARD configuration permit drainage of RO water leaving the RO membrane module, which is conveniently used for laboratory services, such as glass cleaning, bain-marie operation, autoclave refilling. Demineralizers are delivered with a 9 liter pressurized storage tank equipped with a level control to stop water purification when the maximal tank capacity is attained. Larger tanks are available on request.

Optional accessories may be supplied on request, such as a 254 nm UV sterilizer unit or a 0.2 micron filtration unit. Demineralizers in the SIMPLEX configuration feature a two-LED display for conductivity on the control panel.

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Reverse osmosis (RO) is the process that reverses the natural osmotic water flux occurring from a dilute to a concentrated aqueous solution separated by a semi-permeable membrane by the application of a hydrostatic pressure difference across the membrane. By reverse osmosis, the water salt content is reduced on average by 95-99%. Up to 99% of organic matter, bacteria and other impurities is also removed.

The demineralizers of the OSMODEMI series produce high purity water by treatment with reverse osmosis followed by ion-exchange to further reduce water salt content through a bed of ion-exchanging mixed resins. The water produced this way meets the following requirements:

TOC: $< 50 \text{ ppb}$
Water conductivity leaving the RO module: about 2% of that of the feed water (5-10% for the Osmodemi 4-8-12)
Water conductivity leaving the ion-exchange cartridge: $< 0.5 \mu\text{S/cm}$ (i.e., better than $2 \text{ M}\Omega/\text{cm}$).

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The OSMODEMI demineralizers are available either in the SIMPLEX or the STANDARD configuration to meet a broad range of production requirements.

OSMODEMI demineralizers in the SIMPLEX configuration are highly recommended when low/medium amounts of high purity water are required. They are particularly convenient for their

- low capital costs
- simple and easy operation.

Typical applications are in gold or gas production, and in the electronic industry.

OSMODEMI demineralizers in the STANDARD configuration are highly recommended when a broad capacity (from 4 to 60 liter/hour) of RO or high purity water is required. They are particularly convenient for their

- low capital costs
- simple and easy operation
- possibility to drain RO and/or high purity water
- pressurized water storage tanks
- complete instrumentation panel.

Typical applications are in scientific, research and analysis laboratories.

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**Purificatori acqua
ad osmosi inversa e scambio ionico**
**Water purification units
by reverse osmosis and ion-exchange**

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
Equipments


- dual function pre-filter cartridge for dechlorination and 5 micron filtration (code PF 4 – 12)
- reverse osmosis membrane module
- pressurized storage tank with automated stop control driven by water level inside the tank
- de-ionizer cartridge with DF 04 ion-exchange mixed bed resins
- 3-digit electronic digital meter displaying water conductivity leaving the RO membrane module and the ion-exchange cartridge
- drainage stopcocks for RO water and water purified by both RO and ion-exchange
- elegant and sturdy steel cabinet.

Optional accessories:

- UV-ray sterilizer for water post-treatment leaving the ion-exchange cartridge
- 0.2 micron filtration unit
- large storage tanks (e.g., 25 Lt capacity type euro 55; 50 Lt capacity type WM 0120)
- automated water softening unit for feed water pre-treatment.

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
Equipments

- dual function pre-filter cartridge for dechlorination and 5 micron filtration (code PF 4 – 12)
- reverse osmosis membrane module
- de-ionizer cartridge with DF o4 ion-exchange mixed bed resins
- two-LED battery operated conductivity meter
- elegant and sturdy steel cabinet.

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
Equipments

- pre-filtration through an activated charcoal cartridge and a 5 micron filter
- positive-displacement pump
- reverse osmosis membrane module
- de-ionizer cartridge with ion-exchange mixed bed resins
- microprocessor operated electronic controller for safe and fully automatic service. A liquid crystals display shows the different operation phases (in Italian and English)
- automated fluxing at the end of cycle
- elegant and sturdy steel cabinet.

Optional accessories

- storage tank with electric level sensor with the «easy connection» system to facilitate hooking it up to the de-ionizer.

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Equipments

- pre-filtration through an activated charcoal cartridge and a 5 micron filter
- positive-displacement pump
- reverse osmosis membrane module
- flowmeter to measure and show the delivery of the pure water produced
- pressurized storage tank with automated stop control driven by water level inside the tank
- de-ionizer cartridge with ion-exchange mixed bed resins
- electronic controller for safe and fully automatic service.
- 3-digit electronic digital meter displaying water conductivity leaving the RO membrane module and the ion-exchange cartridge
- drainage stopcocks for RO

- water and water purified by both RO and ion-exchange
- possibility to collect water far from the machine
- elegant and sturdy steel cabinet.

Optional accessories

- UV-ray sterilizer for water post-treatment leaving the ion-exchange cartridge
- 0.2 micron filtration unit
- large storage tanks (e.g., 25 Lt capacity type euro 55; 50 Lt capacity type WM 0120)
- automated water softening unit for feed water pre-treatment.

OSMODEMI 20/30 Simplex